

# Disposable Biopsy Forceps

The wide choice of disposable biopsy forceps ensures that you are perfectly equipped for every application. MICRO-TECH offers forceps with diameters of 1.8 mm, 2.3 mm and 3.0 mm in addition to lengths of 120, 180, 230 and 250 cm. Regardless of whether they are tapered, with or without a spike, coated or uncoated and with standard or toothed spoons – all models are characterised by their high reliability. The excellent cutting edge of the biopsy forceps allows you to take diagnostically conclusive tissue samples in a safe, easy manner.

## Performance Characteristics:

- reliability
- very comfortable to use
- diagnostically conclusive biopsies
- wide product variety
- high quality riveted scissors joints
- working channel-friendly design



Toothed spoon



Rivt-hinge design

## Rapid identification:

Jaw section protective cap:

- red – with spike
- yellow – without spike

Coating:

- yellow – paediatric forceps  
Ø 1.8 mm, 180 cm long
- orange – gastric forceps  
Ø 2.3 mm, 180 cm long
- blue – colon forceps  
Ø 2.3 mm, 230 cm long
- orange or blue – macro forceps  
Ø 3.0 mm, 180 und 230 cm long

## Specifications:

REF	Diameter (mm)	Length (cm)	Spike	Covering	Jaw
Forceps with coating					
NBF03-11023180	2.3	180	NO	Yes	standard
NBF03-11023230	2.3	230	NO	Yes	standard
NBF13-11018180	1.8	180	NO	Yes	serrated
NBF13-11023230	2.3	230	NO	Yes	serrated
NBF03-11123180	2.3	180	Yes	Yes	standard
NBF13-11123230	2.3	230	Yes	Yes	serrated