


Fistula-Occluding Esophageal Stents

The fully covered design can occlude the esophageal fistula effectively. And the spherical part in the middle of the stent can prevent side leakage effectively.



Performance Characteristics:

- self-expanding
- Nitinol mesh with atraumatic ends
- excellent positional stability
- high radial force
- resistant and elastic covering
- fully covered stents available
- high radiopacity
- extraction threads for removal and repositioning
- guide wire-compatible up to 0.035 inches

Specifications:

REF	Centre Dia. (mm)	End Dia. (mm)	Length (mm)	Covering (mm)
Stents with end-to-end covering				
ST01-994.14.070	14	20	70	70
ST01-994.16.080	16	22	80	80
ST01-994.18.090	18	24	90	90
ST01-994.20.100	20	26	100	100
ST01-994.22.110	22	28	110	110
ST01-994.24.120	24	30	120	120