


## Pancreatic Pseudocyst Stents



### Performance Characteristics:

- new bridge between gastrointestinal lumen and pancreas
- self-expanding stent with mushroom and umbrella ends, more safety  
(mushroom avoid the stent dropping into the pancreas, umbrella for easy-retrieval)
- special safe lock to ensure the release process in the safe condition
- fully covering, to prevent the hyperplasia
- innovative drainage method, less complication, higher efficiency

## Specifications:

REF	Centre Dia. (mm)	End Dia. (mm)	Covering (mm)	Length (mm)	Specification
ST33-103.10.015	10	20	15	15	Covered mushroom-umbrella
ST33-103.10.020	10	20	20	20	Covered mushroom-umbrella
ST33-103.10.025	10	20	25	25	Covered mushroom-umbrella
ST33-103.10.030	10	20	30	30	Covered mushroom-umbrella
ST33-103.12.015	12	22	15	15	Covered mushroom-umbrella
ST33-103.12.020	12	22	20	20	Covered mushroom-umbrella
ST33-103.12.025	12	22	25	25	Covered mushroom-umbrella
ST33-103.12.030	12	22	30	30	Covered mushroom-umbrella
ST33-103.14.015	14	24	15	15	Covered mushroom-umbrella
ST33-103.14.020	14	24	20	20	Covered mushroom-umbrella
ST33-103.14.025	14	24	25	25	Covered mushroom-umbrella
ST33-103.14.030	14	24	30	30	Covered mushroom-umbrella
ST33-103.16.015	16	26	15	15	Covered mushroom-umbrella
ST33-103.16.020	16	26	20	20	Covered mushroom-umbrella
ST33-103.16.025	16	26	25	25	Covered mushroom-umbrella
ST33-103.16.030	16	26	30	30	Covered mushroom-umbrella