

Hemoclips

The innovative clipping system from MICRO-TECH combines precise positioning with a reliable holding force. This makes the system suitable for effective hemostasis as well as for endoscopic markings. The rotatable clip can be opened and closed multiple times, thus facilitating the exact repositioning of the clip. The clip can be used with any endoscopes equipped with working channels with a minimum diameter of 2.8 mm.

Performance Characteristics:

- Rotatable design for precise placement
- Repeated opening and closing possible
- Excellent holding power
- Wide opening
- Smooth running even under difficult anatomy



Wide range of about 12 mm

Specifications:

REF	opening Diameter (mm)	Working Length (mm)	Sheath Diameter (mm)
ROCC-D-26-195-C	12mm (max)	1950	2.6
ROCC-D-26-230-C	12mm (max)	2300	2.6