

# Injection Needles



Indications for endoscopy to introduce a sclerosing agent or vasoconstrictor into selected sites to control actual or potential bleeding lesions in the digestive system; and the injection of saline to aid in Endoscopic Mucosal Resection (EMR), polypectomy procedures and to control non-variceal hemorrhage.

### Performance Characteristics:

- Lock to ensure that it would not accidentally extend during insertion to damage the scope and to have a stable needle penetration.
- Needle release button design, gently press the button to return the needle fully into the sheath for protecting the endoscope from any damage.
- Ergonomic handle design: complete needle insertion and with drawal with one hand.

### Specifications:

REF	Need Dia. (G)	Needle Length (mm)	Catheter Dia. (mm)	Working Length (mm)
IN12-225231802	22	5	2.3	1800
IN12-225232302	22	5	2.3	2300
IN12-255231802	22	5	2.3	1800
IN12-255232302	22	5	2.3	2300