


# Stone Extraction Baskets



This product is used for endoscopic capture and extraction of biliary calculi.

## Performance Characteristics:

- Specially designed cable to ensure good shape retention.
- Autorotation during basket withdrawal allows for easier capture and better stone retention.

## Specifications:

REF	Basket Diameter (mm)	Sheath Diameter (mm)	Working Length (mm)	Description
SRB-T-9/12-20	20	2.4	2000	4 Wires, Hexagonal
SRB-T-12/15-20	25	2.4	2000	4 Wires, Hexagonal
SRB-T-15/18-20	30	2.4	2000	4 Wires, Hexagonal
SRB-T-18/21-20	35	2.4	2000	4 Wires, Hexagonal