



Biliary Stents



The Biliary Stent is preloaded on the delivery system for ERCP or PTCD. The stent for ERCP is placed through the working channel of the endoscope. The system is introduced into the biliary duct through the papilla under visual control.

Performance Characteristics:

- self-expanding
- available for TTS and PTCD
- Nitinol wire mesh with atraumatic ends
- release under visual endoscopic control
- enormous position stability
- high radiopacity
- guiding wire passable to a maximum of 0.035 inches



Specifications:

REF	Centre Dia. (mm)	End Dia. (mm)	Length (mm)	Covering (mm)
Stents Without Covering				
ST03-001.10.040	10	10	40	Without
ST03-001.10.060	10	10	60	Without
ST03-001.10.080	10	10	80	Without
ST03-001.10.100	10	10	100	Without
Stents With Partial Covering				
ST03-002.10.040	10	10	40	30
ST03-002.10.060	10	10	60	50
ST03-002.10.080	10	10	80	70
ST03-002.10.100	10	10	100	90
Stents With end-to-end Covering				
ST03-003.10.040	10	10	40	40
ST03-003.10.060	10	10	60	60
ST03-003.10.080	10	10	80	80