

# Disposable Biopsy Forceps

The wide choice of disposable biopsy forceps ensures that you are perfectly equipped for every application. The excellent cutting edge of the biopsy forceps allows you to take diagnostically conclusive tissue samples in a safe, easy manner.



## Performance Characteristics:

- reliability
- very comfortable to use
- diagnostically conclusive biopsies
- wide product variety
- high quality riveted scissors joints
- working channel-friendly design

## Specifications:

REF	Diameter (mm)	Working Length (cm)	Spike	Coating	Jaws
Forceps without coating					
NBF11-11018180	1.8	180	no	without	serrated
NBF01-11118120	1.8	120	yes	without	standard
Forceps with coating					
NBF02-11018120	1.8	120	no	grey	standard
NBF02-11118180	1.8	180	yes	yellow	standard