



Y-Shape Tracheal Stents

A world innovation.

MICRO-TECH is the world's first company to present a self-expanding Y-stent for bypassing respiratory tract stenoses in the area around the carina. Its innovative design ensures a stable position and guarantees significant treatment advances in the trachea and in both main bronchi.



Performance Characteristics:

- Y-design
- self-expanding
- Nitinol mesh with atraumatic ends
- excellent positional stability
- high radial force
- resistant and elastic covering
- high radiopacity
- compatible for 2 guide wires of up to 0.035 inches

Specifications:

REF	Trachea Dia. (mm)	Trachea Length (mm)	Main Bronchi Dia. (mm)	Main Bronchi Length (mm) Right/Left	Covering (no covering) Right/Left (mm)
Stent without covering					
ST02-961.20.045	20	45	14	20/30	without
Stent with Partial covering					
ST02-962.16.040	16	40	12	20/30	20(5)/30
ST02-962.20.050	20	50	14	20/30	20(5)/30
Stent with end-to-end covering					
ST02-964.20.050	20	50	14	20/30	20/30