

Disposable Cleaning Brushes

The disposable cleaning brushes from MICRO-TECH ensure a perfectly cleaned endoscope, thus reducing the risk of cross-contamination to a minimum.



Performance Characteristics:

- 2 brush heads for time-saving cleaning
- protective beads on the brush heads
- brushes made of resistant nylon
- individually hygienically packed

Specifications:

REF	Catheter Dia. (mm)	Brushes Dia. (mm)	Working Length (mm)
Cleaning Brushes with one head			
CB18-0050/18	1.8	5	1800
CB18-T0050/23	1.8	5~7 Tapered	2300
Cleaning Brushes with two heads			
CB16-5050/23	1.6	5 / 5	2300
CB18-3030/10	1.8	3 / 3	1000
CB18-3030/18	1.8	3 / 3	1800
CB18-3030/23	1.8	3 / 3	2300
CB18-5050/23	1.8	5 / 5	2300
CB18-6060/23	1.8	6 / 6	2300
CB18-T5050/23	1.8	5~7 / 5~7 Tapered	2300
CB18-T50110/23	1.8	5~7 Tapered / 11	2300
CB23-5050/23	2.3	5 / 5	2300
CB23-T5050/23	2.3	5~7 / 5~7 Tapered	2300
Cleaning Brushes set			
CB18-T0050/23-B	1.8	5~7 Tapered	2300
* CB13-50120/01	1.3	5 / 12	-
CB18-T5050/23-A	1.8	5~7 / 5~7 Tapered	2300
* CB13-00120/00	1.3	12	-
CB18-T5050/23-B	1.8	5~7 / 5~7 Tapered	2300
* CB13-50120/01	1.3	5 / 12	-